

**QUALITY CONTROL CERTIFICATE**

**Cat. No. MAB320P**  
**Dengue NS1 Monoclonal Antibody**

<b>Host</b>	Mouse
<b>Clone Determination</b>	13B06
<b>Immunogen</b>	DENV 2 NS1 Recombinant Protein
<b>Immunoglobulin Class</b>	IgG1-kappa
<b>Specificity</b>	Dengue NS1: serovar 1, 2, 3, & 4 recognition
<b>Purification</b>	Protein A chromatography from <i>in vitro</i> production
<b>Lot Number</b>	
<b>% Purity</b>	> 95% pure determined by size-exclusion chromatography
<b>pI</b>	6.5 – 7.1
<b>Concentration</b>	
<b>Buffer</b>	1 x PBS pH 7.6
<b>Preservatives</b>	None
<b>Expiration Date</b>	2 years from date of manufacture
<b>Quantity</b>	
<b>Volume</b>	
<b>Applications</b>	ELISA pairs with MAB321P, suitable for Lateral Flow
<b>Storage</b>	-20°C

**FOR RESEARCH AND *IN VITRO* USE ONLY**

**MAINE BIOTECHNOLOGY SERVICES, INC.**  
1037R FOREST AVENUE, PORTLAND, ME 04103 USA ~ Tel.: 800-925-9476 / 207-797-5454 / Fax: 207- 797-5595  
Email: sales@mainebiotechnology.com ~ www.mainebiotechnology.com