

Glycyl-tRNA Synthetase (EJ)

Antigen Specification

Product Number: 11100

Description:

Human Glycyl-tRNA Synthetase (EC 6.1.1.14; EJ).
Recombinant antigen for in vitro research and manufacturing use only.

Immunological function:

Binds IgG-type human auto-antibodies.

Origin:

Recombinant. Expressed by recombinant baculovirus (*Autographa californica* multiple nuclear polyhedrosis virus; AcMNPV) infection of *Spodoptera frugiperda* Sf9 insect cells.

Expression construct:

Full-length cDNA coding for the human Glycyl-tRNA Synthetase fused to a deca-histidine purification tag.

Biochemical tests:

SDS-PAGE (purity > 80%); Western blot with i: anti-EJ autoantibody-positive sample; ii: monoclonal anti-His-tag antibody.

Calculated molecular weight:

79 kDa

Calculated isoelectric point:

pH 6.17

Immunological tests/Functionality:

Immunodot analyses with positive/negative samples.

Recommended buffer/storage and handling conditions:

Recommendations for storage buffer: neutral to slightly alkaline pH and 20% glycerol as cryoprotective agent.
Storage conditions: -70°C or below.
Repeated freeze/thaw cycles should be avoided.

Copyright 2013 - 2021 DIARECT GmbH