

<b>PROCESS:</b> Product Specifications	<b>ISSUE DATE:</b> 07.02.2023	<b>DOC NO:</b> 201553BBI CT-PRDS-61
<b>TITLE:</b> Trypsin	<b>REVISION:</b> 11	<b>PAGE:</b> 1 of 1
<b>AUTHORISED:</b> Laboratory Supervisor		<b>APPROVED:</b> SQR Leader

**Product Code:** 201553BBI      **EC no.:** 3.4.21.4      **EDQM Certificate no.:** R1-CEP 2001-209-Rev03

Aspect	Description	Units	Limits	Method
<b>Form</b>	White to off-white, salt-free, freeze-dried powder	N/A	N/A	N/A
<b>Preparation</b>	Bovine pancreas	N/A	N/A	N/A
<b>Activity</b>	Trypsin	U/mg material	≥ 2500	BP6.2.2-WI20
<b>Solubility</b>	Deionised water	50 000 U/ml	Soluble	BP6.2.2-WI105
<b>Additional Data</b>	Loss on Drying (4 hours @ 60 °C)	% m/m	≤5.0	BP6.2.2-WI27
<b>Retest Date</b>	Re-assay	Months	48	N/A
<b>Storage</b>	Temperature	°C	-20	N/A
<b>Shipping</b>	Temperature	°C	Ambient	N/A
<b>Packaging</b>	Dependent on volume / weight	N/A	N/A	N/A
<b>Labels</b>	a) <u>Large Label</u> Name of Product Product Code Lot Number Net Weight Manufacture Date Retest Date Storage Temperature Storage in Transit  b) <u>Small Label</u> Name of Product Product Code Lot No. Net Weight Manufacture Date Retest Date Storage Temperature Storage in Transit	N/A             °C °C             °C °C	N/A             -20 Ambient             -20 Ambient	N/A
<b>Samples</b>	Retention	1 g	6	N/A
	Quality Control	1 g	3	
<b>Special Instructions</b>	None			