

CERTIFICATE OF ANALYSIS

Hyaluronidase

Product Source:	Ovine Testes	Lot no./Batch Number:	SXXXXH39
Product Code:	081239H03	Date of Manufacture:	MM-YYYY
EC No.:	3.2.1.35	Retest Date:	MM-YYYY
Shipping Temperature:	≤ -20 °C	Country of Manufacture:	South Africa
Storage Temperature:	-20 °C		
Description:	White to yellowish, amorphous, freeze-dried powder		

HEALTH AND SAFETY

Application: The product is supplied is for in vitro use only. To be used for research and/or further manufacturing purposes and is not to be used in humans or animals.

PRODUCT INFORMATION

TEST	SPECIFICATION	RESULT
Assay	≥400 u/mg	
Loss on drying	≤5 %	
Identification	Positive	
pH in deionised water (3mg/ml)	3.5-6.5	
Total Aerobic Microbial Count	≤1000 cfu/g	
Total Combined Mould and Yeast Count	≤100 cfu/g	
<i>Salmonella</i>	0 Count/10g	
<i>Pseudomonas aeruginosa</i>	0 Count/g	
<i>E.coli</i>	0 Count/g	
Deionised Water (absence of any extraneous matter)	Soluble (10 mg/ml)	

QUALITY ASSURANCE

I confirm that the above batch has been produced in accordance with specifications and standard operating procedures of BBI Solutions Cape Town.

Name:		Position:	
Signed:		Date:	