

## CERTIFICATE OF ANALYSIS

*Hyaluronidase  $\geq$  1000 u/mg*

|                       |   |                         |              |
|-----------------------|---|-------------------------|--------------|
| Product Source:       | Bovine Testes                           | Lot no./Batch Number:   | SXXXXH46     |
| Product Code:         | 081246BBI                               | Date of Manufacture:    | MM-YYYY      |
| EC No.:               | 3.2.1.35                                | Retest Date:            | MM-YYYY      |
| Shipping Temperature: | 2 – 8 °C                                | Country of Manufacture: | South Africa |
| Storage Temperature:  | -20 °C                                  |                         |              |
| Description:          | White to yellowish, freeze-dried powder |                         |              |

### HEALTH AND SAFETY

**Application:** The product is supplied for in vitro use only. To be used for research and/or further manufacturing purposes and is not to be used in humans or animals.

### PRODUCT INFORMATION

| TEST                   | SPECIFICATION    | RESULT |
|------------------------|------------------|--------|
| Assay                  | $\geq$ 1000 u/mg |        |
| Appearance of solution | $\leq$ Ref III   |        |
| pH                     | 4.5-7.5          |        |

### QUALITY ASSURANCE

I confirm that the above batch has been produced in accordance with specifications and standard operating procedures of BBI Solutions Cape Town.

|         |  |           |  |
|---------|--|-----------|--|
| Name:   |  | Position: |  |
| Signed: |  | Date:     |  |