

Gold Conjugate – Data Sheet

Conjugate:	Goat anti Mouse IgG	40nm Gold
Product Code:	BA.GAM40	
Lot Number:	22110099	

Description of test	Document reference	BBI Specification	Result	Pass/Fail
Optical Density at 520nm at manufacture:	QCP064	10.0 – 10.4	10.1	Pass
Mean Particle Diameter:	QCP007	37.0 – 43.0nm	40.6nm	Pass
Coefficient of Variation of Particle Size:	QCP007	≤ 8%	≤ 8%	Pass
Buffer:	Conjugate is suspended in 2mM Sodium Tetraborate at pH 9.0 containing 0.095% Sodium Azide.			
Date of Manufacture:	14 NOV 2022			

Ref: DS007 Issue 6
 Print Date: 15 Nov. 22
 Effective Date: Refer to Q Pulse

Page 1 of 1

USA: 1-207-797-5454 | International: +44 (0)1495 363 000 | info@bbisolutions.com | www.bbisolutions.com