

Product Analysis

EM Goat anti-Mouse IgG (H&L): 10nm Gold

Product Information

Product Code:	EM.GMHL10
Batch Number:	19060090
Protein Absorption:	Minimum cross reaction to Human, Bovine and Horse serum proteins
Protein Concentration:	14µg/ml
Storage Buffer:	Tris Buffered Saline pH8.2 containing 1% BSA, 15mM NaN ₃ (preservative) and 30% Glycerol
Storage:	Store at 2-8°C

Quality Control

Protein Detected:	<10pg after silver enhancement	Passed
Optical Density:	3.2	Passed
Mean Diameter:	10.0nm	Passed
Coefficient of Variation:	7.6%	Passed
Percentage of Single Particles:	96%	Passed

For in vitro research only. Not for diagnostic use

USA: 1-800-423-8199 | International: +44 (0)2920 767 499 | info@bbisolutions.com | www.bbisolutions.com

Part of **BBI Group**

Registered office: c/o Berry Smith LLP, Haywood House, Dumfries Place, Cardiff, CF10 3GA.
BBI Solutions is the trading name of BBI Solutions OEM Limited. Registered in England and Wales Number 8368483