

Product Analysis

EM Goat anti-Human IgG: 20nm Gold

Product Information

Product Code:	EM.GAH20
Batch Number:	19110217
Protein Absorption:	Not absorbed
Protein Concentration:	13.1µg/ml
Storage Buffer:	Tris Buffered Saline pH8.2 containing 1% BSA, 15mM NaN ₃ (preservative) and 30% Glycerol
Storage:	Store at 2-8°C

Quality Control

Protein Detected:	<1.5ng	Passed
Optical Density:	4.0	Passed
Mean Diameter:	19.7nm	Passed
Coefficient of Variation:	3.9%	Passed
Percentage of Single Particles:	93%	Passed

For in vitro research only. Not for diagnostic use