

Product Analysis

EM Bovine Serum Albumin: 15nm Gold

Product Information

Product Code:	EM.BSA15
Batch Number:	014564
Protein Absorption:	Not Applicable
Storage Buffer:	Tris Buffered Saline pH8.2 containing 1% BSA, 15mM NaN ₃ (preservative) and 30% Glycerol
Storage:	Store at 2-8°C

Quality Control

Protein Detected:	Nothing Detected	Passed
Optical Density:	4.0	Passed
Mean Diameter:	14.9nm	Passed
Coefficient of Variation:	3.9%	Passed
Percentage of Single Particles:	97%	Passed

For in vitro research only. Not for diagnostic use