

Product Analysis

EM Goat anti-Human IgG: 20nm Gold

Product Information

Product Code:	EM.GAH20
Batch Number:	024477
Protein Absorption:	Not absorbed
Protein Concentration:	12.6µg/ml
Storage Buffer:	Tris Buffered Saline pH8.2 containing 1% BSA, 15mM NaN ₃ (preservative) and 30% Glycerol
Storage:	Store at 2-8°C

Quality Control

Protein Detected:	<200pg	Passed
Optical Density:	4.1	Passed
Mean Diameter:	20.9nm	Passed
Coefficient of Variation:	5.2%	Passed
Percentage of Single Particles:	100%	Passed

For in vitro research only. Not for diagnostic use