



The Immunogold Specialist

Product Analysis

EM Rabbit anti-Goat IgG (H&L): 10nm Gold

Product Information

Product Code: EM.RAG 10
Batch Number: 14876
Protein Absorption: Not Absorbed
Protein Concentration: 13.88 µg / ml
Storage Buffer: Tris Buffered Saline pH8.2 containing 1% BSA,
15mM NaN₃ (preservative) and 30% Glycerol
Storage: Store at 4°C

Quality Control

Protein Detected:	7pg after silver enhancement	Passed
Optical Density:	3.2	Passed
Mean Diameter:	9.2nm	Passed
Coefficient of Variation:	5.6%	Passed
Percentage Single Particles:	98%	Passed

QC approved: _____

For in vitro research only. Not for diagnostic use

To reorder or for further technical support:

BBInternational

Golden Gate, Ty Glas Avenue,
Llanishen, Cardiff, CF14 5DX, UK

Tel: 44 (0) 2920 747 232

Fax: 44 (0) 2920 747 242

E-mail: info@bbigold.com

www.bb-international.com

***A complete catalogue of products can be found on our website or CD
Rom***