

HD Grade Data Sheet - 20nm Gold Colloid

Information	
Product	Gold colloid – 20nm
Product Code	HD.GC20.OD5
Batch Number	24110227
Storage Buffer	Suspended in H ₂ O, no preservative Residual chemical left from manufacture
Storage Instructions	Store at 2-8°C – DO NOT FREEZE
Date of Manufacture	29 November 2024 (Use within one month of opening)

Technical Data	
Mean Diameter	19.7nm
Coefficient of Variation	≤8%
Optical Density @ 520nm	5.07
Gold Chloride Concentration	0.05%

USA: 1-800-423-8199 | International: +44 (0)2920 767 499 | info@bbisolutions.com | www.bbisolutions.com