

## CERTIFICATE OF ANALYSIS

*Hyaluronidase >350 U/mg*

Product Source:	Bovine Testes	Lot no./Batch Number:	SXXXXH45
Product Code:	081245BBI	Date of Manufacture:	MM-YYYY
EC No.:	3.2.1.35	Retest Date:	MM-YYYY
Shipping Temperature:	Ambient	Country of Manufacture:	South Africa
Storage Temperature:	2 – 8 °C		
Description:	White to yellowish, freeze-dried powder		

### HEALTH AND SAFETY

**Application:** The product is supplied for in vitro use only. To be used for research and/or further manufacturing purposes and is not to be used in humans or animals.

### PRODUCT INFORMATION

TEST	SPECIFICATION	RESULT
Assay	≥350 u/mg	
Total Aerobic Microbial Count (In house)	Record (cfu/g)	
Total Combined Yeast & Mould Count (In house)	Record (cfu/g)	
0.02M Sodium Phosphate Buffer pH 6.8-7.0	Soluble (5 mg/ml)	
Deionised water	Soluble (5 mg/ml)	

### QUALITY ASSURANCE

I confirm that the above batch has been produced in accordance with specifications and standard operating procedures of BBI Solutions Cape Town.

<b>Name:</b>		<b>Position:</b>	
<b>Signed:</b>		<b>Date:</b>	