

Product:	TRYPSIN INHIBITOR
Product code:	SOY2
CAS number:	9035-81-8
EINECS number:	232-906-9
Source:	Soybean
Form:	A freeze-dried material
Storage conditions:	Store desiccated at -15°C or below. Allow to come to room temperature before opening. Before returning to storage, re-desiccate under vacuum over silica gel for a minimum of four hours. Re-seal before returning to -15°C or below.
Unit definition:	That amount of inhibitor that will inhibit one unit of Trypsin activity (BAEE Units)
Activity:	Between 11,000 and 16,000 BAEE units of inhibition/mg material
Solubility:	Dissolves readily at 5mg/ml in 0.067M Potassium phosphate buffer, pH 7.6 to give a clear solution.