

## EM Grade Data Sheet

| Information          |   |
|----------------------|---|
| Product              | Gold colloid – 40nm   |
| Product Code         | EM. GC40  |
| Batch Number         | 21040114  |
| Storage Buffer       | Suspended in H <sub>2</sub> O, no preservative<br>Residual chemical left from manufacture |
| Storage Instructions | Store at 2-8°C – DO NOT FREEZE  |
| Expiry Date          | April 2022 (Use within one month of opening)  |

| Technical Data              |        |
|-----------------------------|--------|
| Mean Diameter               | 40.5nm |
| Coefficient of Variation    | ≤8%    |
| Optical Density @ 520nm     | 3.00   |
| Gold Chloride Concentration | 0.03%  |

USA: 1-800-423-8199 | International: +44 (0)2920 767 499 | [info@bbisolutions.com](mailto:info@bbisolutions.com) | [www.bbisolutions.com](http://www.bbisolutions.com)

Part of **BBI Group**

Registered office: c/o Berry Smith LLP, Haywood House, Dumfries Place, Cardiff, CF10 3GA.  
BBI Solutions is the trading name of BBI Solutions OEM Limited. Registered in England and Wales Number 8368483