

## Gold Conjugate –Data Sheet

<b>Conjugate:</b>	Monoclonal anti-Biotin:	40nm Gold
<b>Product Code:</b>	BA.MAB40	
<b>Lot Number:</b>	21100036	

Description of test	Document reference	BBI Specification	Result	Pass/Fail
<b>Optical Density at 520nm at manufacture:</b>	QCP064	10.0 – 10.4	10.2	Pass
<b>Mean Particle Diameter:</b>	QCP007	37.0 – 43.0nm	40.4nm	Pass
<b>Coefficient of Variation of Particle Size:</b>	QCP007	≤ 8.0%	≤ 8.0%	Pass
<b>Buffer:</b>	Conjugate is suspended in 2mM sodium Tetraborate at pH 7.2 containing 0.095% Sodium Azide.			
<b>Stability of this product is guaranteed when stored at 2-8°C until:</b>	11 OCT 2023			

Ref: DS055 Issue 5  
 Print Date: 12 Oct. 21  
 Effective Date: Refer to QPulse

Page 1 of 1

USA: 1-207-797-5454 | International: +44 (0)1495 363 000| [info@bbisolutions.com](mailto:info@bbisolutions.com) | [www.bbisolutions.com](http://www.bbisolutions.com)

Part of **BBI Group**

Registered office: c/o Berry Smith LLP, Haywood House, Dumfries Place, Cardiff, CF10 3GA.  
 BBI Solutions is the trading name of BBI Solutions OEM Limited. Registered in England and Wales Number 8368483