

Product Analysis

EM Goat anti-Rabbit IgG: 15nm Gold

Product Information

Product Code:	EM.GAR15
Batch Number:	21110049
Protein Absorption:	Minimum cross reaction to Human serum proteins
Protein Concentration:	8.45µg/ml
Storage Buffer:	Tris Buffered Saline pH8.2 containing 1% BSA, 5mM NaN ₃ (preservative) and 30% Glycerol
Storage:	Store at 2-8°C

Quality Control

Protein Detected:	<1.5ng	Passed
Optical Density:	4.1	Passed
Mean Diameter:	14.3nm	Passed
Coefficient of Variation:	6.9%	Passed
Percentage of Single Particles:	95%	Passed

For in vitro research only. Not for diagnostic use