

Product Analysis

EM Streptavidin: 10nm Gold

Product Information

Product Code:	EM.STP10
Batch Number:	22070092
Protein Absorption:	Not Applicable
Protein Concentration:	11.46µg/ml
Storage Buffer:	Tris Buffered Saline pH8.2 containing 1% BSA, 15mM NaN ₃ (preservative) and 30% Glycerol
Storage:	Store at 2-8°C

Quality Control

Protein Detected:	<100pg after silver enhancement	Passed
Optical Density:	3.1	Passed
Mean Diameter:	9.3nm	Passed
Coefficient of Variation:	8.7%	Passed
Percentage of Single Particles:	96%	Passed

A pdf protocol for the use of Streptavidin conjugates is available on request

For in vitro research only. Not for diagnostic use

USA: 1-800-423-8199 | International: +44 (0)2920 767 499 | info@bbisolutions.com | www.bbisolutions.com

Part of **BBI Group**

Registered office: c/o Berry Smith LLP, Haywood House, Dumfries Place, Cardiff, CF10 3GA.
BBI Solutions is the trading name of BBI Solutions OEM Limited. Registered in England and Wales Number 8368483