

Gold Conjugate –Data Sheet

Conjugate:	Monoclonal anti Fluorescein:	40nm Gold
Product Code:	BA.FITC40	
Lot Number:	23030069	

Description of test	Document reference	BBI Specification	Result	Pass/Fail
Optical Density at 520nm at manufacture:	QCP064	10.0 – 10.4	10.2	Pass
Mean Particle Diameter:	QCP007	37.0 – 43.0nm	40.5nm	Pass
Coefficient of Variation of Particle Size:	QCP007	≤ 8.0%	≤ 8.0%	Pass
Buffer:	Conjugate is suspended in 2mM sodium Tetraborate at pH 9.0 containing 0.095% Sodium Azide.			
Stability of this product is guaranteed when stored at 2 – 8°C until:	13 JUN 2023			

Ref: DS069.2
 Print Date: 15 Mar. 23
 Effective Date: Refer to QPulse

Page 1 of 1

USA: 1-207-797-5454 | International: +44 (0)1495 363 000 | info@bbisolutions.com | www.bbisolutions.com

Part of **BBI Group**

Registered office: c/o Berry Smith LLP, Haywood House, Dumfries Place, Cardiff, CF10 3GA.
 BBI Solutions is the trading name of BBI Solutions OEM Limited. Registered in England and Wales Number 8368483