

## Product Analysis

### EM Goat anti-Rabbit IgG: 5nm Gold

## Product Information

Product Code:	EM.GAR5
Batch Number:	19060055
Protein Absorption:	Minimum cross reaction to Human serum proteins
Protein Concentration:	16.9µg/ml
Storage Buffer:	Tris Buffered Saline pH8.2 containing 1% BSA, 15mM NaN <sub>3</sub> (preservative) and 30% Glycerol
Storage:	Store at 2-8°C

## Quality Control

Protein Detected:	<10pg after silver enhancement	<b>Passed</b>
Optical Density:	3.1	<b>Passed</b>
Mean Diameter:	5.8nm	<b>Passed</b>
Coefficient of Variation:	12.9%	<b>Passed</b>
Percentage of Single Particles:	99%	<b>Passed</b>

**For in vitro research only. Not for diagnostic use**