

Product Analysis

EM Goat anti-Rabbit IgG: 5nm Gold

Product Information

Product Code:	EM.GAR5 S002
Batch Number:	20040068
Protein Absorption:	Minimum cross reaction to Human serum proteins
Protein Concentration:	15.1µg/ml
Storage Buffer:	Tris Buffered Saline pH8.2 containing 1% BSA, 15mM NaN ₃ (preservative) and 30% Glycerol
Storage:	Store at 2-8°C

Quality Control

Protein Detected:	<10pg after silver enhancement	Passed
Optical Density:	3.1	Passed
Mean Diameter:	5.9nm	Passed
Coefficient of Variation:	11.2%	Passed
Percentage of Single Particles:	100%	Passed

For in vitro research only. Not for diagnostic use