

## **Product Analysis**

### **EM Streptavidin: 20nm Gold**

## **Product Information**

Product Code:	EM.STP20
Batch Number:	22100038
Protein Absorption:	Not Applicable
Protein Concentration:	5.11µg/ml
Storage Buffer:	Tris Buffered Saline pH8.2 containing 1% BSA, 15mM NaN <sub>3</sub> (preservative) and 30% Glycerol
Storage:	Store at 2-8°C

## **Quality Control**

Protein Detected:	<400pg	<b>Passed</b>
Optical Density:	4.1	<b>Passed</b>
Mean Diameter:	20.nm	<b>Passed</b>
Coefficient of Variation:	<8%	<b>Passed</b>
Percentage of Single Particles:	99%	<b>Passed</b>

**A pdf protocol for the use of Streptavidin conjugates is available on request**

**For in vitro research only. Not for diagnostic use**

USA: 1-800-423-8199 | International: +44 (0)2920 767 499 | [info@bbisolutions.com](mailto:info@bbisolutions.com) | [www.bbisolutions.com](http://www.bbisolutions.com)

Part of **BBI Group**

Registered office: c/o Berry Smith LLP, Haywood House, Dumfries Place, Cardiff, CF10 3GA.  
BBI Solutions is the trading name of BBI Solutions OEM Limited. Registered in England and Wales Number 8368483