

## **Product Analysis**

### **EM Bovine Serum Albumin: 10nm Gold**

#### **Product Information**

Product Code:	EM.BSA10
Batch Number:	014536
Protein Absorption:	Not Applicable
Storage Buffer:	Tris Buffered Saline pH8.2 containing 1% BSA, 15mM NaN <sub>3</sub> (preservative) and 30% Glycerol
Storage:	Store at 2-8°C

#### **Quality Control**

Protein Detected:	Nothing Detected	<b>Passed</b>
Optical Density:	3.1	<b>Passed</b>
Mean Diameter:	10.0nm	<b>Passed</b>
Coefficient of Variation:	3.7%	<b>Passed</b>
Percentage of Single Particles:	98%	<b>Passed</b>

**For in vitro research only. Not for diagnostic use**