

## **Product Analysis**

### **EM Goat anti-Human IgG: 20nm Gold**

## **Product Information**

Product Code:	EM.GAH20
Batch Number:	019463
Protein Absorption:	Not absorbed
Protein Concentration:	10.2µg/ml
Storage Buffer:	Tris Buffered Saline pH8.2 containing 1% BSA, 15mM NaN <sub>3</sub> (preservative) and 30% Glycerol
Storage:	Store at 2-8°C

## **Quality Control**

Protein Detected:	<400pg	<b>Passed</b>
Optical Density:	4.1	<b>Passed</b>
Mean Diameter:	20.2nm	<b>Passed</b>
Coefficient of Variation:	6.1%	<b>Passed</b>
Percentage of Single Particles:	97%	<b>Passed</b>

**For in vitro research only. Not for diagnostic use**