

Product Analysis

EM Streptavidin: 5nm Gold

Product Information

Product Code:	EM.STP5
Batch Number:	019970
Protein Absorption:	Not Applicable
Protein Concentration:	23.4µg/ml
Storage Buffer:	Tris Buffered Saline pH8.2 containing 1% BSA, 15mM NaN ₃ (preservative) and 30% Glycerol
Storage:	Store at 2-8°C

Quality Control

Protein Detected:	<100pg after silver enhancement	Passed
Optical Density:	3.0	Passed
Mean Diameter:	5.5nm	Passed
Coefficient of Variation:	12.7%	Passed
Percentage of Single Particles:	97%	Passed

A pdf protocol for the use of Streptavidin conjugates is available on request

For in vitro research only. Not for diagnostic use