

Product Analysis

EM Goat anti-Rabbit IgG: 20nm Gold

Product Information

Product Code:	EM.GAR20
Batch Number:	021224
Protein Absorption:	Minimum cross reaction to Human serum proteins
Protein Concentration:	15.2µg/ml
Storage Buffer:	Tris Buffered Saline pH8.2 containing 1% BSA, 15mM NaN ₃ (preservative) and 30% Glycerol
Storage:	Store at 2-8°C

Quality Control

Protein Detected:	<750pg	Passed
Optical Density:	4.1	Passed
Mean Diameter:	19.8nm	Passed
Coefficient of Variation:	5.3%	Passed
Percentage of Single Particles:	99%	Passed

For in vitro research only. Not for diagnostic use