

Product Analysis

EM Protein A: 10nm Gold

Product Information

Product Code:	EM.PAG10
Batch Number:	15249
Protein Absorption:	Not Applicable
Protein Concentration:	4.63µg/ml
Storage Buffer:	Tris Buffered Saline pH8.2 containing 1% BSA, 15mM NaN ₃ (preservative) and 30% Glycerol
Storage:	Store at 2-8°C

Quality Control

Protein Detected:	30pg after silver enhancement	Passed
Optical Density:	3.2	Passed
Mean Diameter:	4.9nm	Passed
Coefficient of Variation:	9.6%	Passed
Percentage of Single Particles:	97%	Passed

For in vitro research only. Not for diagnostic use