



Product Analysis

EM Bovine Serum Albumin: 20nm Gold

Product Information

Product Code:	EM.BSA20
Batch Number:	16518
Protein Absorption:	Not Applicable
Storage Buffer:	Tris Buffered Saline pH8.2 containing 1% BSA, 15mM NaN ₃ (preservative) and 30% Glycerol
Storage:	Store at 4°C

Quality Control

Protein Detected:	Nothing Detected	Passed
Optical Density:	4.0	Passed
Mean Diameter:	20.1nm	Passed
Coefficient of Variation:	6.1%	Passed
Percentage of Single Particles:	98 %	Passed

For in vitro research only. Not for diagnostic use

Re-Ordering and Technical Support

Visit www.bbisolutions.com

Tel: +44 (0) 29 2074 7232 Fax: +44 (0) 29 2074 7242 Email: info@bbisolutions.com

Visit BBI's technical Forum for all the latest industry discussions and support at

www.bbisolutions.com