

The Immunogold Specialist

Product Analysis

EM Bovine Serum Albumin: 40nm Gold

Product Information

Product Code: EM.BSA 40
Batch Number: 10324
Protein Absorption: Not Applicable
Storage Buffer: Tris Buffered Saline pH8.2 containing 1% BSA
15mM NaN₃ (preservative) and 20% Glycerol
Storage: Store at 4°C

Quality Control

Protein Detected:	nothing detected	Passed
Optical Density:	5.0	Passed
Mean Diameter:	40.8nm	Passed
Coefficient of Variation:	5.2%	Passed
Percentage Single Particles:	93%	Passed

QC approved: _____

For in vitro research only. Not for diagnostic use