

The Immunogold Specialist

Product Analysis

EM Streptavidin: 30nm Gold

Product Information

Product Code:	EM.STP30
Batch Number:	16065
Protein Absorption:	Not Applicable
Protein Concentration:	2.75 µg / ml
Storage Buffer:	Tris Buffered Saline pH 8.2 containing 1% BSA 15mM NaN ₃ (preservative) and 30% Glycerol
Storage:	Store at 4°C

Quality Control

Protein Detected:	750pg	Passed
Optical Density:	5.1	Passed
Mean Diameter:	30.0nm	Passed
Coefficient of Variation:	6.7%	Passed
Percentage Single Particles:	98%	Passed

QC approved: _____

A pdf protocol for the use of Streptavidin conjugates is available on request.

For in vitro research only. Not for diagnostic use