

**Gold Conjugate –Data Sheet**

<b>Conjugate:</b>	Monoclonal anti Horseradish Peroxidase:	40nm Gold
<b>Product Code:</b>	BA.HRP40	
<b>Lot Number:</b>	24010104	

Description of test	Document reference	BBI Specification	Result	Pass/Fail
<b>Optical Density at 520nm at manufacture:</b>	QCP064	10.0 – 10.4	10.3	Pass
<b>Mean Particle Diameter:</b>	SPEC-007	37.0 – 43.0nm	39.8nm	Pass
<b>Coefficient of Variation of Particle Size:</b>	SPEC-007	≤ 8.0%	≤ 8.0%	Pass
<b>Buffer:</b>	Conjugate is suspended in 2mM sodium Tetraborate at pH 9.0 containing 0.095% Sodium Azide.			
<b>Stability of this product is guaranteed when stored at 2 – 8°C until:</b>	15 APR 2024			

Ref: DS070.3  
 Print Date: 19 Jan. 24  
 Effective Date: Refer to QPulse