

Gold Conjugate –Data Sheet

Conjugate:	Goat anti Human IgG (H+L):	40nm Gold
Product Code:	BA.GAHL40	
Lot Number:	24010183	

Description of test	Document reference	BBI Specification	Result	Pass/Fail
Optical Density at 520nm at manufacture:	QCP064	10.0 – 10.4	10.3	Pass
Mean Particle Diameter:	SPEC-007	37.0 – 43.0nm	39.8nm	Pass
Coefficient of Variation of Particle Size:	SPEC-007	≤ 8.0%	≤ 8.0%	Pass
Buffer:	Conjugate is suspended in 2mM sodium Tetraborate at pH 7.2 containing 0.095% Sodium Azide.			
Stability of this product is guaranteed when stored at 2-8°C until:	18 JAN 2024			

Ref: DS018 Issue 10
 Print Date: 25 Jan. 24
 Effective Date: Refer to QPulse