

## Product Analysis

### EM Goat anti-Mouse IgG: 10nm Gold

## Product Information

Product Code:	EM.GAM10
Batch Number:	24030221
Protein Absorption:	Minimum cross reaction to Human, Bovine and Horse serum proteins
Protein Concentration:	15µg/ml
Storage Buffer:	Tris Buffered Saline pH8.2 containing 1% BSA, 15mM NaN <sub>3</sub> (preservative) and 30% Glycerol
Storage:	Short Term Storage at 2-8°C (<6months) Long Term Storage -20°C (>6months)

## Quality Control

Protein Detected:	<7pg after silver enhancement	<b>Passed</b>
Optical Density:	3.1	<b>Passed</b>
Mean Diameter:	10.1nm	<b>Passed</b>
Coefficient of Variation:	6.5%	<b>Passed</b>
Percentage of Single Particles:	92%	<b>Passed</b>

**For in vitro research only. Not for diagnostic use**

USA: 1-800-423-8199 | International: +44 (0)2920 767 499 | [info@bbisolutions.com](mailto:info@bbisolutions.com) | [www.bbisolutions.com](http://www.bbisolutions.com)

Part of **BBI Group**

Registered office: c/o Berry Smith LLP, Haywood House, Dumfries Place, Cardiff, CF10 3GA.  
BBI Solutions is the trading name of BBI Solutions OEM Limited. Registered in England and Wales Number 8368483