

Product Analysis

EM Goat anti-Rabbit IgG: 5nm Gold

Product Information

Product Code:	EM.GAR5
Batch Number:	24070110
Protein Absorption:	Minimum cross reaction to Human serum proteins
Protein Concentration:	21µg/ml
Storage Buffer:	Tris Buffered Saline pH8.2 containing 1% BSA, 15mM NaN ₃ (preservative) and 30% Glycerol
Storage:	Short Term Storage at 2-8°C (<6months) Long Term Storage -20°C (>6months)

Quality Control

Protein Detected:	<10pg after silver enhancement	Passed
Optical Density:	3.1	Passed
Mean Diameter:	5.2nm	Passed
Coefficient of Variation:	<15%	Passed
Percentage of Single Particles:	100%	Passed

For in vitro research only. Not for diagnostic use