

Product Analysis

EM Streptavidin: 20nm Gold

Product Information

Product Code:	EM.STP20
Batch Number:	22100038
Protein Absorption:	Not Applicable
Protein Concentration:	5.11µg/ml
Storage Buffer:	Tris Buffered Saline pH8.2 containing 1% BSA, 15mM NaN ₃ (preservative) and 30% Glycerol
Storage:	Store at 2-8°C

Quality Control

Protein Detected:	400pg/µL	Passed
Optical Density:	4.1	Passed
Mean Diameter:	20.3nm	Passed
Coefficient of Variation:	<8%	Passed
Percentage of Single Particles:	99%	Passed

A pdf protocol for the use of Streptavidin conjugates is available on request

For in vitro research only. Not for diagnostic use

USA: 1-800-423-8199 | International: +44 (0)2920 767 499 | info@bbisolutions.com | www.bbisolutions.com