

Product Analysis

EM Streptavidin: 20nm Gold

Product Information

Product Code:	EM.STP20
Batch Number:	08128_001
Protein Absorption:	Not Applicable
Protein Concentration:	5.1µg/ml
Storage Buffer:	Tris Buffered Saline pH8.2 containing 1% BSA, 15mM NaN ₃ (preservative) and 30% Glycerol
Storage:	Store at 4°C

Quality Control

Protein Detected:	<750pg	Passed
Optical Density:	4.0	Passed
Mean Diameter:	20.3nm	Passed
Coefficient of Variation:	6.3%	Passed
Percentage of Single Particles:	97%	Passed

A pdf protocol for the use of Streptavidin conjugates is available on request

For in vitro research only. Not for diagnostic use